

Rat PEPCK (Phosphoenolpyruvate Carboxykinase) Sandwich ELISA Kit (STJ150377)

STJ150377

GENERAL INFORMATION


Product Type	Kits
Description	This Rat PEPCK (Phosphoenolpyruvate Carboxykinase) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection Limit	0.156--10ng/mL
Sensitivity	10ng/mL
Detection Limit	0.156--10ng/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	PCK1
Symbol	
Uniprot ID	PCKGC_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081