

Rat SLC (Secondary Lymphoid Tissue Chemokine) Sandwich ELISA Kit (STJ150381)

STJ150381

GENERAL INFORMATION


Product Type	Kits
Short Description	This Rat SLC (Secondary Lymphoid Tissue Chemokine) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Rat

PRODUCT PROPERTIES

Detection	62.5--4000pg/mL
Sensitivity	4000pg/mL
Detection Limit	62.5--4000pg/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	
Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081