

## Human RNASE3 Sandwich ELISA Kit (STJ150447)

STJ150447

### GENERAL INFORMATION


|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Kits   |
| <b>Short Description</b> | This Human RNASE3 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA  |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                        |                              |
|------------------------|------------------------------|
| <b>Detection Limit</b> | 39.063--2500pg/mL            |
| <b>Sensitivity</b>     | 2500pg/mL                    |
| <b>Detection Limit</b> | 39.063--2500pg/mL            |
| <b>Assay Time</b>      |                              |
| <b>Storage</b>         | Store at 4°C for 6 months.   |
| <b>Instruction</b>     |                              |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO). |

### TARGET INFORMATION

|                     |                                  |
|---------------------|----------------------------------|
| <b>Sample Types</b> | Serum, Plasma, Biological Fluids |
| <b>Gene</b>         | <a href="#">RNASE3</a>           |
| <b>Symbol</b>       |                                  |
| <b>Uniprot ID</b>   | <a href="#">ECP_HUMAN</a>        |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081