

Human Ubiquitin Carboxyl Terminal Hydrolase L1 Sandwich ELISA kit (STJ150471)

STJ150471

GENERAL INFORMATION


Product Type	Kits
Short Description	This Human Ubiquitin Carboxyl Terminal Hydrolase L1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Limit	0.781--50ng/mL
Sensitivity	50ng/mL
Detection Limit	0.781--50ng/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	ELANE
Symbol	
Uniprot ID	ELNE_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081