

Mouse High Mobility Group Protein B1 Sandwich ELISA kit (STJ150484)

STJ150484

GENERAL INFORMATION


Product Type	Kits
Short Description	This Mouse High Mobility Group Protein B1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	0.313--20ng/mL
Sensitivity	20ng/mL
Detection Limit	0.313--20ng/mL
Assay Time	
Storage	Store at 4°C for 6 months.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene	CORT
Symbol	
Uniprot ID	CORT_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081