

Mouse NT-proBNP (N-terminal pro-Brain Natriuretic Peptide) Sandwich ELISA Kit (STJ150488)

STJ150488

GENERAL INFORMATION


Product Type	Kits
Short Description	This Mouse NT-proBNP (N-terminal pro-Brain Natriuretic Peptide) Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection Limit	15.625--1000pg/mL
Sensitivity	1000pg/mL
Detection Limit	15.625--1000pg/mL
Assay Time	
Storage Instruction	Store at 4°C for 6 months.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Biological Fluids
Gene Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081