

Anti-CD274 antibody [IHC411] (STJ160016)

STJ160016

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal anti-CD274 antibody is suitable for use in Immunohistochemistry.
Applications	IHC-p
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IHC411
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	1:100-1:200
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	
Storage Instruction	Store at 2-8°C. Do not freeze.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
Uniprot ID	PD1L1_HUMAN
Immunogen	Human PD-L1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunohistochemical stain on human lung.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081