

Anti-IgG1 antibody (PE) (STJ16100001)

STJ16100001

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IgG1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	NA
Uniprot ID	IGG1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	