

Anti-HLA-b antibody [JOAN-1] (STJ16100009)

STJ16100009

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-b is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC, IHC-F
Host/Source	Mouse
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	JOAN-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3106
Gene Symbol	HLA-B
Uniprot ID	HLAB_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

