

Anti-CD11c antibody [FK24] (STJ16100012)

STJ16100012

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11c is suitable for use in Immunofluorescence and Immunohistochemistry research applications.
Applications	IF, IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FK24
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IF (1:5) IHC (1:5)
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3687
Gene Symbol	ITGAX
Uniprot ID	ITAX_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	