

Anti-CR3 antibody [BEAR-1] (STJ16100018)

STJ16100018

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11b/CR3 is suitable for use in Immunohistochemistry, Flow Cytometry and Immunoprecipitation research applications.
Applications	IHC-F, FC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BEAR-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3684
Gene Symbol	ITGAM
Uniprot ID	ITAM_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	