

Anti-CR3 antibody [BEAR-1] {FITC} (STJ16100019)

STJ16100019

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11b/CR3 is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F, FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BEAR-1
Concentration	
Conjugation	FITC
Purification	
Dilution Range	FC. 5 μ g is the recommended amount of reagent required for 1×10^6 cells.
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3684
Gene Symbol	ITGAM
Uniprot ID	ITAM_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	