

Anti-CD11a antibody [87-6F9] (STJ16100022)

STJ16100022

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11a is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC, IHC-F, IF, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	87-6F9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 0.5-1.0 Mu g/million cells in 0 1 ml IHC – Frozen 1.0-2.0 Mu g/ml IF 1.0-2.0 Mu g/ml Functional Assays 0 02-2 0 Mu g/ml without azide
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3683
Gene Symbol	ITGAL
Uniprot ID	ITAL_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	