

## Anti-EGFR antibody [2, 00E+09] (STJ16100028)

STJ16100028

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Epidermal Growth Factor Receptor/EGFR is suitable for use in Immunofluorescence, Immunoprecipitation and EM research applications.
<b>Applications</b>	IF, IP, EM
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	2, 00E+09
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IF (1:25-1:100) IP EM
<b>Formulation</b>	PBS with 1% BSA & 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">1956</a>
<b>Gene Symbol</b>	<a href="#">EGFR</a>
<b>Uniprot ID</b>	<a href="#">EGFR_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	