

Anti-CD23 antibody [Tù 1] {FITC} (STJ16100033)

STJ16100033

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody Anti-CD23 is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Tù 1
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	0.5 ml of affinity purified antibody, labelled with FITC, in PBS containing 0.1% sodium Azide and 4mg/ml BSA.
Isotype	IgG3
Storage Instruction	2-8°C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	