

Anti-CD98 antibody [MEM-108] (STJ16100036)

STJ16100036

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD98 is suitable for use in Flow Cytometry, Immunoprecipitation and Immunohistochemistry research applications.
Applications	FC, IP, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-108
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 4 µg/ml-IP to be determined-IHC (paraffin) 10 µg/ml
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence