

## Anti-CD81 antibody [1.3.3.22] (STJ16100041)

STJ16100041

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD81 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
<b>Applications</b>	FC, IHC-F, IHC-P, IF, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1.3.3.22
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC 1-2 µg/million cells in 0 1 ml at 4C-IF 1-2 µg/ml-WB 1-2 µg/ml-IHC (paraffin) Recom
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">975</a>
<b>Gene Symbol</b>	<a href="#">CD81</a>
<b>Uniprot ID</b>	<a href="#">CD81_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	