

## Anti-CD84 antibody [152-1D5] (STJ16100042)

STJ16100042

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD84 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	FC, IHC-F, IF, IP, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	152-1D5
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC 1-2 µg/million cells in 0 1 ml at 4C-IF 1-2 µg/ml-IP 1-2 µg per 100-500 µg of t
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">8832</a>
<b>Gene Symbol</b>	<a href="#">CD84</a>
<b>Uniprot ID</b>	<a href="#">SLAF5_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	