

Anti-CD29 antibody [Moon-4] (STJ16100043)

STJ16100043

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD29 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Moon-4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 0 5-1 0 µg/million cells in 0 1 ml-IF 0 5-1 0 µg/ml-IHC (frozen) 1-2 µg/ml for 30-60 min at RT
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence