

Anti-CD28 antibody [204.12] (STJ16100044)

STJ16100044

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD28 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human, Mouse, Pig, Sheep, Cow

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	204.12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 0 5-1 0 µg/million cells in 0 1 ml-IF 0 5-1 0 µg/ml-IHC (frozen) 1-2 µg/ml for 30-60 min at RT
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	940
Gene Symbol	CD28
Uniprot ID	CD28_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	