

Anti-HLA-DR antigen antibody [LN-3] (STJ16100045)

STJ16100045

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-DR antigen is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LN-3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	serumfree culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence