

Anti-HLA/DR antibody [HL-39] (STJ16100046)

STJ16100046

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA/DR is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HL-39
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 2 µg/ml-IP to be determined
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	