

Anti-CD44 antibody [MEM-85] {PE} (STJ16100064)

STJ16100064

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD44 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-85
Concentration	
Conjugation	PE
Purification	
Dilution Range	FC 1:50
Formulation	PBS with sucrose, BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
Uniprot ID	CD44_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	