

Anti-CD14 antibody [Tük 4] (STJ16100069)

STJ16100069

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody Anti-CD14 is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Tük 4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.5 ml purified antibody, in Phosphate buffered saline solution containing 0.1% sodium Azide.
Isotype	IgG2a
Storage Instruction	2-8°C

TARGET INFORMATION

Gene ID	929
Gene Symbol	CD14
Uniprot ID	CD14_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	