

Anti-CD45RB antibody [MEM-55] {FITC} (STJ16100087)
STJ16100087

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD45RB is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Primates, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-55
Concentration	
Conjugation	FITC
Purification	
Dilution Range	FC using 20 μ l reagent/100 μ l of whole blood or 10 ⁶ cells in a suspension.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081