

## Anti-CD55 antibody [MEM-118] (STJ16100088)

STJ16100088

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD55/DAF is suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	FC, IP, IHC-P, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Primates

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MEM-118
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Tris buffered saline (TBS) solution with 15 mM sodium azide
<b>Isotype</b>	IgM
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	1604
<b>Gene Symbol</b>	CD55
<b>Uniprot ID</b>	DAF_HUMAN
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	