

Anti-CD59 antibody [MEM-43/5] (STJ16100089)

STJ16100089

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD59 is suitable for use in Western Blot, Flow Cytometry, Immunoprecipitation and Immunohistochemistry research applications.
Applications	WB, FC, IP, IHC-P
Host/Source	Mouse
Reactivity	Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-43/5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	966
Gene Symbol	CD59
Uniprot ID	CD59_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081