

Anti-CD47 antibody [BRIC126] (STJ16100098)

STJ16100098

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD47 is suitable for use in Western Blot, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	WB, FC, IHC-F, IP
Host/Source	Mouse
Reactivity	Human, Pig, Sheep, Dog, Bovine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRIC126
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:10
Formulation	TRIS buffered glycine
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
Uniprot ID	CD47_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	