

Anti-CD51 antibody [23C6] (STJ16100100)

STJ16100100

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD51 is suitable for use in Immunofluorescence, Immunoprecipitation, Flow Cytometry and Immunohistochemistry research applications.
Applications	IF, IP, FC, IHC-F
Host/Source	Mouse
Reactivity	Human, Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	23C6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:40-1:80
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence