

Anti-CD15 antibody [MEM-158] (STJ16100105)

STJ16100105

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD15 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-158
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	