

Anti-CD13 antibody [Tük 1] {R-PE} (STJ16100107)

STJ16100107

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody Anti-CD13 is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Tük 1
Concentration	
Conjugation	R-PE
Purification	
Dilution Range	
Formulation	0.5 ml purified antibody, conjugated to R-PE, in Phosphate Buffer Saline) solution containing sucrose, 4 mg/ml BSA and 0.1% sodium azide.
Isotype	IgG1
Storage Instruction	2-8°C

TARGET INFORMATION

Gene ID	290
Gene Symbol	ANPEP
Uniprot ID	AMPN_HUMAN
Immunogen	Human CD13
Immunogen Region	
Specificity	
Immunogen Sequence	