

Anti-CD16 antibody [3G8] {PE} (STJ16100112)

STJ16100112

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD16 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3G8
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 4mg/ml BSA and 0.1% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	