

Anti-CD19 antibody [SJ25-C1] {FITC} (STJ16100114)

STJ16100114

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD19 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SJ25-C1
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 4mg/ml BSA and 0.1% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	930
Gene Symbol	CD19
Uniprot ID	CD19_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	