

Anti-PECAM-1 antibody [BV8] (STJ16100122)

STJ16100122

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-PECAM-1/CD31 is suitable for use in Flow Cytometry, ELISA, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, ELISA, IF, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BV8
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	5175
Gene Symbol	PECAM1
Uniprot ID	PECA1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	