

## Anti-HLA-Class I antibody [W6/32] (STJ16100127)

STJ16100127

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-HLA-Class I is suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot, Immunohistochemistry, Immunocytochemistry and ELISA research applications.
<b>Applications</b>	FC, IP, WB, IHC-F, ICC, ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Cow, Primates, Human, Cat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	W6/32
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution</b>	
<b>Range</b>	
<b>Formulation</b>	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	