

Anti-HLA-Class I antibody [W6/32] {Biotin} (STJ16100128)

STJ16100128

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-Class I is suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot, Immunohistochemistry, Immunocytochemistry and ELISA research applications.
Applications	FC, IP, WB, IHC-F, ICC, ELISA
Host/Source	Mouse
Reactivity	Cow, Primates, Human, Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	W6/32
Concentration	
Conjugation	Biotin
Purification	
Dilution	
Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	