

**Anti-HLA-Class I antibody [W6/32] {PE} (STJ16100130)**  
STJ16100130**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-HLA-Class I is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Cow, Primates, Human, Cat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	W6/32
<b>Concentration</b>	
<b>Conjugation</b>	PE
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
<b>Isotype</b>	IgG2a
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen Sequence</b>