

Anti-Campylobacter jejuni antibody [57-24] (STJ16100136)

STJ16100136

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Campylobacter jejuni is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Bacteria

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	57-24
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide and 0.1% BSA
Isotype	IgG3
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	