

Anti-SLPI antibody [31] (STJ16100138)

STJ16100138

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SLPI is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, ELISA, IP, IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	31
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50 IP 1:50 WB 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6590
Gene Symbol	SLPI
Uniprot ID	SLPI_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	