

Anti-CD4-antigen antibody [15-8A2] (STJ16100150)

STJ16100150

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD4-antigen is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F, FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	15-8A2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:100
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence