

Anti-CD8 antibody [15-11C5] (STJ16100151)

STJ16100151

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD8 is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F, FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	15-11C5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:100
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	