

Anti-B3 Integrin/CD61 antibody [BV4] (STJ16100154)

STJ16100154

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B3 Integrin/CD61 is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, ELISA, IP, IHC-P, WB
Host/Source	Mouse
Reactivity	Human, Bovine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BV4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 ELISA 1:50 WB 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	