

Anti-CD1b antibody [100-1A5] (STJ16100159)

STJ16100159

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD1b is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	ELISA, FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	100-1A5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	910
Gene Symbol	CD1B
Uniprot ID	CD1B_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	