

Anti-blood group antigen A1B antibody [HE-24] (STJ16100161)

STJ16100161

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-blood group antigen A1B is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HE-24
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence

Clone	Agglutination of tested blood groups				
	A ₁	B	A ₁ B	A ₂ B	O
HE-193 (anti-A)	+	-	+	+	-
HE-195 (anti-A)	+	-	+	+	-
HEB-29 (anti-B)	-	+	+	+	-
HE-24 (anti-A ₁ B)	+	-	+	-	-

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081