

Anti-Myeloid Specific Antigen antibody [BM-3] (STJ16100170)

STJ16100170

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Myeloid Specific antigen is suitable for use in Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	IHC-F, IF, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BM-3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC frozen sections (1-2 µg/m) IF (0.5-1.0 µg/ml) and IP.
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

