

Anti-Myeloid Specific Antigen antibody [BM-4] (STJ16100171)

STJ16100171

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Myeloid Specific antigen is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-F, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BM-4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC frozen sections (1-2 µg/m) IF (0.5-1.0 µg/ml).
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

