

## Anti-TLR9 antibody [5G5] {Biotin} (STJ16100184)

STJ16100184

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-TLR9 is suitable for use in research applications.
<b>Applications</b>	
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Universal

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5G5
<b>Concentration</b>	
<b>Conjugation</b>	Biotin
<b>Purification</b>	
<b>Dilution Range</b>	FC 1:50 ELISA 1:50 WB 1:50 IP 1:50
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">54106</a>
<b>Gene Symbol</b>	<a href="#">TLR9</a>
<b>Uniprot ID</b>	<a href="#">TLR9_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	