

Anti-JAM-1 antibody [BV12] (STJ16100185)

STJ16100185

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-JAM-1 is suitable for use in Flow Cytometry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	FC, IF, IP, WB
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BV12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
Uniprot ID	JAM1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	