

Anti-CD105 antibody [MEM-226] (STJ16100186)

STJ16100186

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD105 is suitable for use in Flow Cytometry, Immunoprecipitation and Western Blot research applications.
Applications	FC, IP, WB
Host/Source	Mouse
Reactivity	Rat, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-226
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	