

Anti-CD105 antibody [MEM-226] {FITC} (STJ16100187)

STJ16100187

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD105 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-226
Concentration	
Conjugation	FITC
Purification	
Dilution Range	FC (20 µl reagent/100 µl sample).
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence